


TriCollect

The clever solution that avoids plastic waste

Sustainable sample transportation from sender to laboratory



Reduce your costs

Eliminate costs of plastic bags

Avoid disposal of plastic waste as potentially infectious material

Free your staff from unpacking plastic bags

Ensure secure processes

Work safely, in accordance with the packing instruction P 650 (UN 3373)

Ensure safe handling and transportation of samples and request forms due to coherent and patented concept

Support complete coagulation of serum tubes

Improve your workflow

Act environmentally friendly without effort

Enjoy a well-structured process from tube collection to arrival in the laboratory

Reduce time for unpacking of samples and request forms

